


Powerful 3D printer Planmeca Creo™



Brand new!



From vision to reality

Planmeca Creo™ is our new and powerful 3D printer for creating a wide range of indications. It is the perfect tool for dental labs and large clinics aiming to expand their production capabilities and increase efficiency.

The 3D printer uses Digital Light Processing (DLP) technology to create objects out of UV curable resin. It has its own dedicated software, which is part of all deliveries.

Perfect
your craft

Create your vision
with Planmeca Creo™



Planmeca Creo™

Technical specifications

Print technology	Digital Light Processing (DLP)
Resolution (XY)	30–75 µm, 1200 x 1920
Layer thickness	5–200 µm
Light source	LED
Guaranteed LED life	10,000 hours
Build volume	130 x 90 x 130 mm (5.1 x 3.5 x 5.1 in.)
AC input	120/240 V
Weight	27 kg (59.5 lbs)
Dimensions (W x H x D)	42 x 32 x 72 cm (16.5 x 12.5 x 28 in.)
Minimum required clearances	Front: 42 cm (16.5 in.) Sides: 10 cm (3.9 in.) Rear: 5 cm (2 in.) Top: 10 cm (3.9 in.)
Operating temperature	15–32 °C (59–89.6 °F)
Materials from	Planmeca

Indications

Dental models	YES
Dental splints	YES
Surgical guides	YES
Casts	YES
Trays	YES
Temporary fillings	Coming soon
Ortho models	Coming soon
Bases	Coming soon
Gingiva masks	Coming soon

PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, DentoVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumo, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca mRomexis, Planmeca Noma, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanMill, Planmeca Planosil, Planmeca PlanPure, Planmeca PlanScan, Planmeca ProCeph, Planmeca ProFace, Planmeca ProLED, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca Serenus, Planmeca SingLED, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Verity, Planmeca Waterline Cleaning System, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartPan, SmartTouch, Trendy and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

